




Naturest® Kenkopad®

Transform your regular mattress into a uniquely designed, adaptive sleep environment

- Non-allergenic latex foundation moulds to your body's contours for better comfort and support
- Cocooning revolutionary RAM™ (radial-axis magnetism) Technology
- Made from ecologically sound and renewable non-allergenic latex
- Gradient-density nodules provide a gentle massage-like stimulation
- Latex foundation wicks away moisture and promotes improved ventilation
- Natural materials meet fire regulations without flame retardants or other chemicals



Advanced relaxation technology

The Nikken Naturest® Kenkopad® is designed from latex foam sculpted into an 'egg-crate' surface with specially treated tips. These tips provide a relaxing massage effect while sleeping and allow for air movement and dissipation of the body's natural perspiration during the night. The latex also yields to pressure for an ideal degree of sleep comfort and is naturally hygienic due to its resistance to microbial growth and dust mites.

Using Nikken's unique RAM™ technology magnetic arrays, the Naturest® Kenkopad® is complemented with a specially designed magnetic field for relaxation and improved sleep quality.

Natural, renewable materials

The Naturest® Kenkopad® is made from blended non-allergenic latex foam. The cover is woven from all-natural, fully renewable cotton, wool,

and bamboo fibres, produced without chemical pesticides or insecticides. These natural materials improve the breathability, assist with the dissipation of moisture, and consequently help create a drier, more comfortable, sleeping environment.

How to use

Place the Naturest® Kenkopad® on top of your regular mattress with the nodules facing upwards. Make sure it is placed so that the tag and zipper are at the top where your head will rest.

Care instructions

Spot clean only, using mild detergent and warm water. Do not oversaturate. Air-dry at room temperature away from sunlight or excessive heat. If folding the Naturest Kenkopad is necessary, fold with the flat side outward. Do not attempt to fold with the nodules facing out.

| | |
|------------------------|--|
| Technology | RAM™ Magnetic Technology |
| Magnet strength | 3,200 Gauss |
| Materials | Synthetic latex: 38%; natural latex: 34%; graphite: 24%; wool batting: 4% |
| Magnetic Discs | 13/19/26 per pad depending on size, providing 3,200 gauss strength per Naturest® Kenkopad® |
| Outer cover | Natural fibres; a mix of organic cotton and bamboo. |

1110 Naturest® Kenkopad® 74x185

1111 Naturest® Kenkopad® 90x190

1112 Naturest® Kenkopad® 135x190

1113 Naturest® Kenkopad® 150x200